



Cd81 Antibody, HRP conjugated

Product Code	CSB-PA004960DB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35762
Immunogen	Recombinant Mouse CD81 antigen protein (116-201AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD81 antigen (26 kDa cell surface protein TAPA-1) (Target of the antiproliferative 1) (CD antigen CD81), Cd81, Tapa1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd81