



serB Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA018938HB01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0AGB0 |
| Immunogen | Recombinant Escherichia coli Phosphoserine phosphatase protein (1-322AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Phosphoserine phosphatase (PSP) (PSPase) (EC 3.1.3.3) (O-phosphoserine phosphohydrolase), serB |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | serB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |