



serB Antibody, Biotin conjugated

Product Code	CSB-PA018938HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AGB0
Immunogen	Recombinant Escherichia coli Phosphoserine phosphatase protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoserine phosphatase (PSP) (PSPase) (EC 3.1.3.3) (O-phosphoserine phosphohydrolase), serB
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	serB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.