



SCP Antibody, Biotin conjugated

Product Code	CSB-PA321237HD01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14348
Immunogen	Recombinant Epstein-Barr virus Small capsomere-interacting protein (31-176AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small capsomere-interacting protein, SCP
Immunogen Species	Epstein-Barr virus
Research Area	Others
Target Names	SCP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.