



IptA Antibody, FITC conjugated

Product Code	CSB-PA360010HC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ADV1
Immunogen	Recombinant Escherichia coli Lipopolysaccharide export system protein LptA protein (28-185AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipopolysaccharide export system protein LptA, IptA, yhbN
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	IptA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.