



Adrb2 Antibody

Product Code	CSB-PA001392LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10608
Immunogen	Recombinant Rat Beta-2 adrenergic receptor protein (330-418AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-2 adrenergic receptor (Beta-2 adrenoreceptor) (Beta-2 adrenoceptor), Adrb2, Adrb2r
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Adrb2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.