



Reg3b Antibody

Product Code	CSB-PA326510LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25031
Immunogen	Recombinant Rat Regenerating islet-derived protein 3-beta protein (27-175AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regenerating islet-derived protein 3-beta (REG-3-beta) (Pancreatitis-associated protein 1) (Peptide 23) (REG-2) (Regenerating islet-derived protein III-beta) (Reg III-beta) [Cleaved into: Regenerating islet-derived protein 3-beta 16.5 kDa form; Regenerating islet-derived protein 3-beta 15 kDa form], Reg3b, Pap Pap1 Reg2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Reg3b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.