



# Reg3b Antibody

<b>Product Code</b>	CSB-PA326510LA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P25031
<b>Immunogen</b>	Recombinant Rat Regenerating islet-derived protein 3-beta protein (27-175AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Regenerating islet-derived protein 3-beta (REG-3-beta) (Pancreatitis-associated protein 1) (Peptide 23) (REG-2) (Regenerating islet-derived protein III-beta) (Reg III-beta) [Cleaved into: Regenerating islet-derived protein 3-beta 16.5 kDa form; Regenerating islet-derived protein 3-beta 15 kDa form], Reg3b, Pap Pap1 Reg2
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Others
<b>Target Names</b>	Reg3b