




GSTK1 Antibody

Product Code	CSB-PA009978HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Q3
Immunogen	Recombinant Human Glutathione S-transferase kappa 1 protein (2-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase kappa 1 (EC 2.5.1.18) (GST 13-13) (GST class-kappa) (GSTK1-1) (hGSTK1) (Glutathione S-transferase subunit 13), GSTK1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GSTK1

Image



Western blot

All lanes: GSTK1 antibody at 2µg/ml

Lane 1: 293T whole cell lysate

Lane 2: Hela whole cell lysate

Lane 3: A431 whole cell lysate

Lane 4: COLO205 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 26, 32, 25, 21 kDa

Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.