



EMC9 Antibody, FITC conjugated

Product Code	CSB-PA897482LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3B6
Immunogen	Recombinant Human ER membrane protein complex subunit 9 protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ER membrane protein complex subunit 9 (Protein FAM158A), EMC9, C14orf122, FAM158A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EMC9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.