



plp Antibody

Product Code	CSB-PA018202LA01OEI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79826
Immunogen	Recombinant Oncorhynchus mykiss Myelin proteolipid protein (150-218AA)
Raised In	Rabbit
Species Reactivity	Oncorhynchus mykiss
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin proteolipid protein (PLP) (DM20) (Lipophilin), plp
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Research Area	Others
Target Names	plp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.