





AGRP Antibody

Product Code	CSB-PA05699A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00253
Immunogen	Recombinant Human Agouti-related protein (83-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Agouti-related protein, AGRP, AGRT ART
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AGRP

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA05699A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA05699A0Rb at dilution of 1:100