





GADD45G Antibody

Product Code	CSB-PA009164HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95257
Immunogen	Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 gamma protein (1-159AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth arrest and DNA damage-inducible protein GADD45 gamma (Cytokine-responsive protein CR6) (DNA damage-inducible transcript 2 protein) (DDIT-2), GADD45G, CR6 DDIT2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GADD45G

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA009164HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA009164HA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.