



IFNL1 Antibody, Biotin conjugated

Product Code	CSB-PA13249D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU54
Immunogen	Recombinant Human Interferon lambda-1 protein (21-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon lambda-1 (IFN-lambda-1) (Cytokine Zcyto21) (Interleukin-29) (IL-29), IFNL1, IL29 ZCYTO21
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNL1