



QPCTL Antibody, FITC conjugated

Product Code	CSB-PA868351LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXS2
Immunogen	Recombinant Human Glutaminyl-peptide cyclotransferase-like protein (212-382AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaminyl-peptide cyclotransferase-like protein (EC 2.3.2.5) (Golgi-resident glutaminyl-peptide cyclotransferase) (isoQC) (gQC), QPCTL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	QPCTL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.