





ALKBH3 Antibody

Product Code	CSB-PA847209HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q83
Immunogen	Recombinant Human Alpha-ketoglutarate-dependent dioxygenase alkB homolog 3 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-ketoglutarate-dependent dioxygenase alkB homolog 3 (EC 1.14.11.54) (Alkylated DNA repair protein alkB homolog 3) (hABH3) (DEPC-1) (Prostate cancer antigen 1), ALKBH3, ABH3 DEPC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ALKBH3


Image



Immunohistochemistry analysis of human testis tissue using CSB-PA847209HA01HU at dilution of 1:100



Immunocytochemistry analysis of human skeletal muscle tissue using CSB-PA847209HA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA847209HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.