



Tp53 Antibody

Product Code	CSB-PA024077LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02340
Immunogen	Recombinant Mouse Cellular tumor antigen p53 protein (1-387AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellular tumor antigen p53 (Tumor suppressor p53), Tp53, P53 Trp53
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Tp53