



Tp53 Antibody, HRP conjugated

Product Code	CSB-PA024077LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02340
Immunogen	Recombinant Mouse Cellular tumor antigen p53 protein (1-387AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellular tumor antigen p53 (Tumor suppressor p53), Tp53, P53 Trp53
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Tp53
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.