



Tp53 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA024077LC01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02340 |
| Immunogen | Recombinant Mouse Cellular tumor antigen p53 protein (1-387AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cellular tumor antigen p53 (Tumor suppressor p53), Tp53, P53 Trp53 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Cell Biology |
| Target Names | Tp53 |