




Mug1 Antibody

Product Code	CSB-PA328826LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28665
Immunogen	Recombinant Mouse Murinoglobulin-1 protein (700-910AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Mug1

Image



Western Blot

Positive WB detected in: Mouse Liver tissue lysate

All lanes: Mug1 antibody at 1:500

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 166 kDa

Observed band size: 166 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.