



Major pollen allergen Car b 1 isoforms 1A and 1B Antibody, HRP conjugated

Product Code	CSB-PA330956LB01CEC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38949
Immunogen	Recombinant Carpinus betulus Major pollen allergen Car b 1 isoforms 1A and 1B protein (2-160AA)
Raised In	Rabbit
Species Reactivity	Carpinus betulus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Car b 1 isoforms 1A and 1B (Allergen Car b I) (allergen Car b 1)
Immunogen Species	Carpinus betulus
Research Area	Others