



IFNT1 Antibody, FITC conjugated

Product Code	CSB-PA322795LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15696
Immunogen	Recombinant Bovine Interferon tau-1 protein (24-195AA)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon tau-1 (IFN-tau-1) (Antiluteolysin) (Trophoblast antiluteolytic protein) (Trophoblast protein 1) (TP-1) (Trophoblastin), IFNT1
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	IFNT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.