



nrdB Antibody, Biotin conjugated

Product Code	CSB-PA303206HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69924
Immunogen	Recombinant Escherichia coli Ribonucleoside-diphosphate reductase 1 subunit beta protein (2-376AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonucleoside-diphosphate reductase 1 subunit beta (EC 1.17.4.1) (Protein B2) (Protein R2) (Ribonucleotide reductase 1), nrdB, ftsB
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	nrdB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.