



ecoRVM Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA366147LA01ENL |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04393 |
| Immunogen | Recombinant Escherichia coli Modification methylase EcoRV protein (1-298AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Modification methylase EcoRV (M.EcoRV) (EC 2.1.1.72) (Adenine-specific methyltransferase EcoRV), ecoRVM |
| Immunogen Species | Escherichia coli |
| Research Area | Others |
| Target Names | ecoRVM |