



ecoRVM Antibody

Product Code	CSB-PA366147LA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04393
Immunogen	Recombinant Escherichia coli Modification methylase EcoRV protein (1-298AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase EcoRV (M.EcoRV) (EC 2.1.1.72) (Adenine-specific methyltransferase EcoRV), ecoRVM
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	ecoRVM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.