



mpt51 Antibody, HRP conjugated

Product Code	CSB-PA14937B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WQN7
Immunogen	Recombinant Mycobacterium tuberculosis MPT51/MPB51 antigen protein (27-299AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MPT51/MPB51 antigen, mpt51, fbpC1 fbpD mpb51
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	mpt51
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.