




APRT Antibody

Product Code	CSB-PA001954DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07741
Immunogen	Recombinant Human Adenine phosphoribosyltransferase protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenine phosphoribosyltransferase (APRT) (EC 2.4.2.7), APRT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	APRT

Image


Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA001954DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.