




BCL7C Antibody

Product Code	CSB-PA845147DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUZ0
Immunogen	Recombinant Human B-cell CLL/lymphoma 7 protein family member C protein (58-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	B-cell CLL/lymphoma 7 protein family member C, BCL7C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCL7C

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA845147DSR2HU at dilution of 1:100