





BOC Antibody

Product Code	CSB-PA866299ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWV1
Immunogen	Recombinant Human Brother of CDO protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Brother of CDO (Protein BOC), BOC
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	BOC

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA866299ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA866299ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.