




DHRS2 Antibody

Product Code	CSB-PA615674ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13268
Immunogen	Recombinant Human Dehydrogenase/reductase SDR family member 2, mitochondrial protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dehydrogenase/reductase SDR family member 2, mitochondrial (EC 1.1.1.-) (Dicarbonyl reductase HEP27) (Protein D) (Short chain dehydrogenase/reductase family 25C member 1), DHRS2, SDR25C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DHRS2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA615674ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.