





DPP8 Antibody

Product Code	CSB-PA744278ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6V1X1
Immunogen	Recombinant Human Dipeptidyl peptidase 8 protein (10-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl peptidase 8 (DP8) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 1) (DPRP-1) (Dipeptidyl peptidase VIII) (DPP VIII) (Prolyl dipeptidase DPP8), DPP8, DPRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPP8

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA744278ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA744278ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.