





GSTO2 Antibody

Product Code	CSB-PA872498ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4Y5
Immunogen	Recombinant Human Glutathione S-transferase omega-2 protein (124-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase omega-2 (GSTO-2) (EC 2.5.1.18) (Glutathione S-transferase omega 2-2) (GSTO 2-2) (Glutathione-dependent dehydroascorbate reductase) (EC 1.8.5.1) (Monomethylarsonic acid reductase) (MMA(V) reductase) (EC 1.20.4.2), GSTO2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSTO2

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA872498ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA872498ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.