





IGFBP4 Antibody

Product Code	CSB-PA011098ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22692
Immunogen	Recombinant Human Insulin-like growth factor-binding protein 4 protein (22-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-like growth factor-binding protein 4 (IBP-4) (IGF-binding protein 4) (IGFBP-4), IGFBP4, IBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IGFBP4

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA011098ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA011098ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.