




IL22 Antibody

Product Code	CSB-PA867113ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZX6
Immunogen	Recombinant Human Interleukin-22 protein (34-179AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-22 (IL-22) (Cytokine Zcyto18) (IL-10-related T-cell-derived-inducible factor) (IL-TIF), IL22, ILTIF ZCYTO18
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL22

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA867113ESR1HU at dilution of 1:100