




# ANGPTL4 Antibody

<b>Product Code</b>	CSB-PA866314ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BY76
<b>Immunogen</b>	Recombinant Human Angiotensin-related protein 4 protein (26-406AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Angiotensin-related protein 4 (Angiotensin-like protein 4) (Hepatic fibrinogen/angiotensin-related protein) (HFARP), ANGPTL4, ARP4 HFARP PGAR
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	ANGPTL4

**Image**


Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA866314ESR1HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.