



ITGA2B Antibody, Biotin conjugated

Product Code	CSB-PA011865LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08514
Immunogen	Recombinant Human Integrin alpha-IIb protein (639-887AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-IIb (GPalph IIb) (GPIIb) (Platelet membrane glycoprotein IIb) (CD antigen CD41) [Cleaved into: Integrin alpha-IIb heavy chain; Integrin alpha-IIb light chain, form 1; Integrin alpha-IIb light chain, form 2], ITGA2B, GP2B ITGAB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ITGA2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.