




H2-D1 Antibody

Product Code	CSB-PA365684LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01900
Immunogen	Recombinant Mouse H-2 class I histocompatibility antigen, D-D alpha chain protein (25-311AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	H-2 class I histocompatibility antigen, D-D alpha chain (H-2D(D)), H2-D1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	H2-D1

Image



Western blot
 All lanes: H2-D1 antibody at 2µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 42 kDa
 Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.