





HYAL2 Antibody

Product Code	CSB-PA618635ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12891
Immunogen	Recombinant Human Hyaluronidase-2 protein (367-473AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronidase-2 (Hyal-2) (EC 3.2.1.35) (Hyaluronoglucosaminidase-2) (Lung carcinoma protein 2) (LuCa-2), HYAL2, LUCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HYAL2

Image



Western blot
 All lanes: Hyaluronidase-2 antibody at 4µg/ml +
 Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 54 kDa
 Observed band size: 54 kDa



Immunohistochemistry of paraffin-embedded
 human kidney tissue using CSB-
 PA618635ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.