




IL13RA2 Antibody

Product Code	CSB-PA618901ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14627
Immunogen	Recombinant Human Interleukin-13 receptor subunit alpha-2 protein (27-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-13 receptor subunit alpha-2 (IL-13 receptor subunit alpha-2) (IL-13R subunit alpha-2) (IL-13R-alpha-2) (IL-13RA2) (Interleukin-13-binding protein) (CD antigen CD213a2), IL13RA2, IL13R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL13RA2

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA618901ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.