





COPS2 Antibody

Product Code	CSB-PA005789ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61201
Immunogen	Recombinant Human COP9 signalosome complex subunit 2 protein (1-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	COP9 signalosome complex subunit 2 (SGN2) (Signalosome subunit 2) (Alien homolog) (JAB1-containing signalosome subunit 2) (Thyroid receptor-interacting protein 15) (TR-interacting protein 15) (TRIP-15), COPS2, CSN2 TRIP15
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	COPS2

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA005789ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005789ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.