





KIR2DS4 Antibody

Product Code	CSB-PA012362ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43632
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS4 protein (22-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DS4 (CD158 antigen-like family member I) (MHC class I NK cell receptor) (Natural killer-associated transcript 8) (NKAT-8) (P58 natural killer cell receptor clones CL-39/CL-17) (p58 NK receptor CL-39/CL-17) (CD antigen CD158i), KIR2DS4, CD158I KKA3 NKAT8
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DS4

Image



Western blot
 All lanes: Killer cell immunoglobulin-like receptor 2DS4 antibody at 4μg/ml
 Lane 1: K562 whole cell lysate
 Lane 2: Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 34 kDa
 Observed band size: 34 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA012362ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.