




NEK2 Antibody

Product Code	CSB-PA015699ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51955
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek2 protein (1-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek2 (EC 2.7.11.1) (HSPK 21) (Never in mitosis A-related kinase 2) (NimA-related protein kinase 2) (NimA-like protein kinase 1), NEK2, NEK2A NLK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NEK2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA015699ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.