





TNFSF18 Antibody

Product Code	CSB-PA891791ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNG2
Immunogen	Recombinant Human Tumor necrosis factor ligand superfamily member 18 protein (72-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 18 (Activation-inducible TNF-related ligand) (AITRL) (Glucocorticoid-induced TNF-related ligand) (hGITRL), TNFSF18, AITRL GITRL TL6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFSF18

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA891791ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA891791ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.