




CDC5L Antibody

Product Code	CSB-PA858712ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99459
Immunogen	Recombinant Human Cell division cycle 5-like protein (563-802AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle 5-like protein (Cdc5-like protein) (Pombe cdc5-related protein), CDC5L, KIAA0432 PCDC5RP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CDC5L

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA858712ESR2HU at dilution of 1:100