



CPM Antibody

Product Code	CSB-PA005897ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14384
Immunogen	Recombinant Human Carboxypeptidase M protein (204-443AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase M (CPM) (EC 3.4.17.12), CPM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CPM