



LAT Antibody

Product Code	CSB-PA012767ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43561
Immunogen	Recombinant Human Linker for activation of T-cells family member 1 protein (28-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Linker for activation of T-cells family member 1 (36 kDa phospho-tyrosine adapter protein) (pp36) (p36-38), LAT
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LAT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.