



IRAK4 Antibody

Product Code	CSB-PA011812ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWZ3
Immunogen	Recombinant Human Interleukin-1 receptor-associated kinase 4 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 receptor-associated kinase 4 (IRAK-4) (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-64), IRAK4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRAK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.