



# Ctsg Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA006190LC01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P28293
<b>Immunogen</b>	Recombinant Mouse Cathepsin G protein (21-261AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Relevance</b>	This vimentin-specific protease may regulate the reorganization of vimentin filaments, occurring during cell differentiation, movement and mitosis.
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cathepsin G (EC 3.4.21.20) (Vimentin-specific protease) (VSP), Ctsg
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	Ctsg