



Ctsg Antibody, Biotin conjugated

Product Code	CSB-PA006190LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28293
Immunogen	Recombinant Mouse Cathepsin G protein (21-261AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Relevance	This vimentin-specific protease may regulate the reorganization of vimentin filaments, occurring during cell differentiation, movement and mitosis.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin G (EC 3.4.21.20) (Vimentin-specific protease) (VSP), Ctsg
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Ctsg