




MORF4L1 Antibody

Product Code	CSB-PA883368ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU8
Immunogen	Recombinant Human Mortality factor 4-like protein 1 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mortality factor 4-like protein 1 (MORF-related gene 15 protein) (Protein MSL3-1) (Transcription factor-like protein MRG15), MORF4L1, MRG15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MORF4L1

Image



Western blot

All lanes: Mortality factor 4-like protein 1


antibody at 2µg/ml + A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 42, 38, 27 kDa

Observed band size: 42 kDa



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA883368ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.